

BVI

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¹ Ophthalmologica Review: Retinal Toxicity of Medical Devices Used during Vitreoretinal Surgery: A Critical Overview acute retinal damage after using a toxic perfluorocarbon for vitreoretinal surgery- retina 37:1140-1151, 2017.

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IOBA | Cell Culture Lab: Preclinical In Vitro
Cytotoxicity study of PFO medical devices (-ARCADOPHTA XII-)

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Study Report

In vitro cytotoxicity study of PFO medical devices on cultures of human retinal pigment epithelial cells and porcine neuroretina -ARCADOPHTA-

Sponsor: **ARCADOPHTA** **11 Rue Antoine Ricord**
31100 Toulouse
France

Conclusions:

The product "PFO medical device" lot numbers 310832 and 400929 (A and B) were provided by the company ARCADOPHTA and tested with RPE cell cultures,

- Lots 310832 and 400929 are NON-CYTOTOXIC.**

It was concluded based on the differences in terms of viability compared to the controls used in the experiments performed following the UNE-EN ISO 10993-5:2009 "Tests for in vitro cytotoxicity".

The product "PFO medical device" lot numbers 310832 and 400929 (A and B) were provided by the company ARCADOPHTA and tested with porcine neuroretina cultures,

- Lots 310832 and 400929 are NON-CYTOTOXIC.**

It was concluded based on the differences in terms of the cell/tissue degeneration compared to the controls used in the experiments performed following the UNE-EN ISO 10993-5:2009 "Tests for in vitro cytotoxicity".

Version: **01** **Date:** **26/10/2020**

Study Director: **JC Pastor Jimeno**

Experiments Directors: **Cristina Andrés Iglesias**
Iván Fernandez Bueno

Facilities Director: **MJ Maldonado López**

Quality Assurance Unit: **MA García Hontoria**

This report shall not be reproduced except in full without the specific written approval of the IOBA, Valladolid, Spain
- The test results relate only to the items tested -
IOBA - Campus Miguel Delibes - Paseo de Belén 17 - 47011 - Valladolid
Tel. 983 42 35 59 | Fax 983 18 37 23 | ioba@ioba.med.uva.es | www.ioba.es

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